

KT
MEDICAL



KT Guide Wire

廣泰導引線
衛部醫器製字第006144號

KT Medical

廣泰導引線 (康聚醫學科技股份有限公司(廣泰集團))

KT Medical Inc.(Kuang Tai Group)

We developed the very advanced technology to manufacture our guidewire series to meet doctor's expectation and adopted the highest standard of quality control system to make sure quality sustained. With the most cutting-edge technology we developed our own PTFE coating process that ensures the coating is applied strongly and smoothly on the product which ensures our product has the lowest risk of being damaged during operations. The smoothly coated surface of our PTFE product helps the doctor to operate our KT wires to move and turn as precise as a well-trained sniper, never miss a target. Our GW are designed to meet the most common requests from the market while we are able to custom make preferred properties such as tip shape and stiffness of the wire tip.


我們的導線產品是利用自行開發出來的最先進技術製程以達到醫師在手術上的需求，並且以高規格要求品質水準能夠維持一致。利用最先進的PTFE塗層技術可以確保塗層能夠有足夠的包覆強度以及滑順的表面，讓醫生在使用的時候降低導線受損的機率，並且能夠順暢精確地操作。目前公司產品以市場主要規格為主，但我們仍能提供客製化規格(導線頭端形狀以及軟硬度) 以符合客戶的需求。

Wire Tip 導線頭端

With J-type or Straight type tip, doctors can easily and safely turn the tip to find its path in the lengthy procedure. The rounded tip-end also provides safety in the tortuous anatomy.

利用J型頭端或者直線型頭端，醫生可以輕易以及安全地操作手術，而圓順的頂端可以降低手術時傷害血管的風險。








Stiffness 軟硬度

KT Wire is designed with 3 Levels of stiffness: Soft, standard, and stiff for doctor's preference which are achieved by different core wire design. This ensures doctors get the most suitable material and have the confidence using the wire for the procedure.

KT導線利用心線設計來達到不同柔軟度，分別為Stiff、Standard 以及 Soft，可依照醫師的喜好自由選擇。



PTFE Coating - pre-coated wires PTFE預塗層技術

Our PTFE coils are made with pre-coated wires. The process ensures the integrity of the coating coverage and better adhesion. This advantage eliminates the concern of coating damages.

我們的PTFE線圈是利用預塗層的線材所繞製而成，這種製程的好處在於表面塗層的完整性以及包覆性優於傳統噴塗工法，降低操作過程中塗層剝落的風險。



● Single End KT Guidewire-J Tip 廣泰導引線 - 單邊柔軟J頭

OD (Inch)	L(cm)	Tip Description	Flexibility	Order No.
0.035	80	3.0mm J tip	Standard	PSJS-35-080-3-std
0.035	90	3.0mm J tip	Standard	PSJS-35-090-3-std
0.035	150	3.0mm J tip	Soft	PSJS-35-150-3-soft
0.035	150	3.0mm J tip	Standard	PSJS-35-150-3-std
0.035	150	6.0mm J tip	Standard	PSJS-35-150-6-std
0.035	150	3.0mm J tip	Stiff	PSJS-35-150-3-stiff
0.035	180	3.0mm J tip	Standard	PSJS-35-180-3-std
0.035	260	3.0mm J tip	Standard	PSJS-35-260-3-std
0.038	150	3.0mm J tip	Soft	PSJS-38-150-3-soft
0.038	150	3.0mm J tip	Standard	PSJS-38-150-3-std
0.038	150	3.0mm J tip	Stiff	PSJS-38-150-3-stiff
0.038	180	3.0mm J tip	Standard	PSJS-38-180-3-std
0.038	260	3.0mm J tip	Standard	PSJS-38-260-3-std

● Single End KT Guidewire-Straight Tip 廣泰導引線 - 單邊柔軟直頭

OD (Inch)	L(cm)	Tip Description	Flexibility	Order No.
0.035	80	Straight tip	Stiff	PSIS-35-080-stiff
0.035	90	Straight tip	Stiff	PSIS-35-090-stiff
0.035	150	Straight tip	Soft	PSIS-35-150-soft
0.035	150	Straight tip	Standard	PSIS-35-150-std
0.035	150	Straight tip	Stiff	PSIS-35-150-stiff
0.038	150	Straight tip	Standard	PSIS-38-150-std
0.038	150	Straight tip	Stiff	PSIS-38-150-stiff
0.035	180	Straight tip	Standard	PSIS-35-180-std
0.035	180	Straight tip	Stiff	PSIS-35-180-stiff
0.038	180	Straight tip	Stiff	PSIS-38-180-stiff
0.035	260	Straight tip	Standard	PSIS-35-260-std
0.035	260	Straight tip	Stiff	PSIS-35-260-stiff
0.038	260	Straight tip	Stiff	PSIS-38-260-stiff

● Double End KT Guidewire-J/Straight Tip 廣泰導引線 - 雙邊柔軟J/直頭

OD (Inch)	L(cm)	Tip Description	Flexibility	Order No.
0.035	80	3.0 J/Straight	Soft	PSJF-35-080-3-soft
0.035	80	3.0 J/Straight	Standard	PSJF-35-080-3-std
0.035	90	3.0 J/Straight	Soft	PSJF-35-090-3-soft
0.035	90	3.0 J/Straight	Standard	PSJF-35-090-3-std
0.035	150	3.0 J/Straight	Standard	PSJF-35-150-3-std
0.035	260	3.0 J/Straight	Standard	PSJF-35-260-3-std
0.038	150	3.0 J/Straight	Stiff	PSJF-38-150-3-stiff
0.038	260	3.0 J/Straight	Standard	PSJF-38-260-3-std

All J tip are finger-straightenable. More options are available. Please contact us for more details.

所有J頭端皆可單手操作。若您需求之規格不在上述表格內，歡迎與我們聯絡。

業務及技術諮詢

For any inquires, please contact:

康聚醫學科技股份有限公司
KT Medical Inc.
T : 886-7-6955155 ext. 7612 or 7623
F : 886-7-6955166
E : medsales@kuangtai.com

製造廠址

Manufacturing Site:

康聚醫學科技股份有限公司
KT Medical Inc.
高雄市路竹區路科一路9號
No. 9, Luke 1st Rd., Luzhu Dist.,
Kaohsiung City 82151, Taiwan



MD 732304